

## How to Use and Search "The LIST" & Safety Topics

Clicking "The LIST & Articles" on the main menu of [mrisafety.com](http://mrisafety.com) brings you to a page like this:

The screenshot shows the MRISafety.com search results page. At the top left is the logo 'MRISafety.com' with a red 'A' above it. Below the logo is a navigation bar with links: 'Research:' (B), 'The LIST' (C), 'Safety Topics' (D), and 'Home' (E). To the right is a search field (F) with a magnifying glass icon and the text '(whole or partial words)'. Below the search field is a pagination bar (G) with page numbers 1 through 10, 'Next', and 'Last'. Below the pagination bar is a display count: 'Displaying 1 - 20 of 4844' with a dropdown menu set to '20'. The main content is a table (I) with three columns: 'Object Description' (with an upward arrow), 'Object Status', and 'Safety Topic / Subject'. The first row (H) has a magnifying glass icon (H) and a red 'J' next to the description. The 'Object Status' column for the first two rows is highlighted in yellow (K) and contains 'Conditional 5 More...' and 'Conditional 7 More...' respectively. The last two rows have 'Unsafe 1 More...' in the 'Object Status' column. The table contains the following data:

Object Description ↑	Object Status	Safety Topic / Subject
2D Helical, 35 Fibered Platinum Coil. Boston Scientific, <a href="http://www.bostonscientific.com">www.bostonscientific.com</a>	1.5, 3 Conditional 5 More...	Coils, Filters, Stents, and G
3/4" Socket Wrench 3/4" 41-mm Newmatic Medical, <a href="http://www.newmaticmedical.com">www.newmaticmedical.com</a>	3 Conditional 7 More...	Miscellaneous Implants and
30 Caliber, 762 x 39, Copper Jacketed Round, Armor Piercing, Norinco	1.5 Unsafe 1 More...	Miscellaneous Implants and
357 Magnum Revolver Model 66-3 Misc. Smith and Wesson	1.5 Unsafe 1 More...	Miscellaneous Implants and

We have labeled important areas of this page as A-K, and will describe each area below:

- A. Is the header of the page. It is a link ([MRISafety.com](http://MRISafety.com)) that will bring you back to the Homepage.
- B. Research is a link that will bring you to the page you see above.
- C. Menu choice: "The LIST" which also brings you to the page you see above.
- D. Menu choice: "Safety Topics" which displays a searchable list of MRI related safety articles.
- E. Menu choice: "Home" takes you to the [MRISafety.com](http://MRISafety.com) homepage.  
**? help** - Clicking "help" displays this Help PDF file in another browser window.
- F. Search field. Here you can type anything and, after press ENTER or clicking the magnifying glass, records that contain what you typed will show up in the list. You can type one word or multiple words - whole words or partial words. More on this function later. Note: search is not case sensitive.
- G. Pages of data available based upon your search criteria. You can move quickly through many pages of data by clicking a page number.
- H. Click this magnifying glass to view the full Object record. (although, nearly all info about Objects is already being displayed on the list screen)
- I. Clicking on column titles will sort the list (toggle either ascending or descending) by the data in that column.
- J. Clicking on the Object description will also take you to the full page Object description.
- K. Clicking "More..." will always show you **more** related information.

Below is an example of searching The LIST (Safety Topic searching works the same way):

The screenshot shows the top navigation bar of 'The LIST' website. It includes links for 'The LIST', 'Safety Topics', and 'Home'. A search bar contains the text 'accu shun' and has a magnifying glass icon and a clear button (X). Below the search bar, it says 'Displaying 1 - 6 of 6' and a dropdown menu is set to '20'. The main content area is a table with three columns: 'Object Description', 'Object Status', and 'Safety Topic / S'. The table contains three rows of search results, each with a magnifying glass icon in the first column.

Object Description ↑		Object Status	Safety Topic / S
<a href="#">AccuDrain External CSF Drainage System(s)</a> Cerebrospinal Fluid (CSF) <b>Shunt</b> Valves and Accessories Integra NeuroSciences www.integralife.com	3	Safe More...	Cerebrospinal F
<a href="#">Accura II</a> Adult Integral <b>Shunt</b> System Phoenix Biomedical Corporation Valley Forge, PA	1.5	Safe More...	Cerebrospinal F
<a href="#">Accura II</a> Pediatric Integral <b>Shunt</b> System Phoenix Biomedical Corporation Valley Forge, PA	1.5	Safe More...	Cerebrospinal F

**Using Search:** You can type one or more words in the Search box - and these can be whole words or partial words. As you can see above, the word fragments "accu shun" have been entered in the search field. This search displays six records that contain BOTH of those word fragments. Clicking "X" on the right side of the search box clears the box and resets the list, making it ready to accept a new search.

Clicking on the "Safety Topics" menu item brings you to a screen like this:

MRIsafety.com

Research: [The LIST](#) [Safety Topics](#) [Home](#)

1 2 3 4 5 6 7 8

Displaying 1 - 20 of 162

	<input type="checkbox"/>	Objects Toggle	Safety Topic / Subject
<input type="checkbox"/>	<input type="checkbox"/>	41	Breast Tissue Expanders and Implants
<input type="checkbox"/>	<input type="checkbox"/>	0	Guidelines to Prevent Excessive Heating and Burns Associated with MRI
<input type="checkbox"/>	<input type="checkbox"/>	61	Cardiac Pacemakers, Implantable Cardioverter Defibrillators (ICDs), and Cardiac
<input type="checkbox"/>	<input type="checkbox"/>	93	Cardiovascular Catheters, Other Catheters, and Accessories
<input type="checkbox"/>	<input type="checkbox"/>	5	Carotid Artery Vascular Clamps
<input type="checkbox"/>	<input type="checkbox"/>	3	Cochlear Implants

The Safety Topics list contains an extra **little pink toggle button on each row**. When you click it, the list expands to display all of the Objects that are related to this Topic/Subject. Below you can see there are three Objects in The LIST that are related to Cochlear Implants. **Click the little pink button again to hide the list of related Objects (it's a toggle switch).**

<input type="checkbox"/>	<input type="checkbox"/>	93	Cardiovascular Catheters, Other Catheters, and Accessories
<input type="checkbox"/>	<input type="checkbox"/>	5	Carotid Artery Vascular Clamps
<input type="checkbox"/>	<input type="checkbox"/>	3	Cochlear Implants

THE List [The LIST](#)

Displaying 1 - 3 of 3

Object Description ↑	Object Status
<input type="checkbox"/> Alpha (M) Magnetic Implant, Bone Conduction Hearing System Sophono, Inc., www.sophono.com	1.5,3 Conditional 5 More...
<input type="checkbox"/> Esteem hearing implant Envoy Medical www.envoymedical.com	1.5 Unsafe 1 More...
<input type="checkbox"/> Nucleus Freedom, Cochlear Implant Cochlear Corporation, www.cochlear.com	1.5 Conditional 5 More...

<input type="checkbox"/>	<input type="checkbox"/>	840	Coils, Filters, Stents, and Grafts
--------------------------	--------------------------	-----	------------------------------------